

MEDIC BAG



E-mail: contact@resqsupply.com Phone: +45 911 19 911



MEDIC BACKPACK



HEMOSTATIC MODULE



INFUSION MODULE



VENTILATION MODULE



INSTRUMENT MODULE



BANDAGING MODULE



STRETCHER MODULE





NO	MODULE	NAME	SPECFICATIONS	QTY	PIECE	FUNCTION	PICTURE
1		Hemostatic sponge module	40 × 30 × 8mm 1 pcs/ bagΦ80 × 60× 8mm 1 pcs/ bagΦ18×8mm 5 pcs/bag	3	bu	Used for emergency situations to control severe surface bleeding and achieve hemostasis.	
2		Wound sponge dressing	120× 80 × 8mm 1 pcs/ bag	3	рс	For temporary external use to control deeper and lage surface bleeding.	
3		Spinning tourniquet	3.8cm × 95cm	2	рс	One hand application. Occlude blood flow of an extremity in the event of a traumatic wound with significant hemorrhage.	
4		TPE tourniquet	2.5cm × 45cm	3	rl	Disposable tear resistant tourniquets provide a secure and comfortable hold on any sized limb.	0
5		Emergency hemostatic bandage	10cm × 45cm	2	rl	Suitable for acute superficial wounds on the body, providing emergency hemostasis, as well as dressing and fixation.	



NO	MODULE	NAME	SPECFICATIONS	QTY	PIECE	FUNCTION	PICTURE
6		Emergency bandage	150mm × 1200mm	2	rl	Multi-functional bandage to treat a wide range of injuries in a compact package.	
7		Compressed gauze	11cm × 350cm	2	рс	An ultra-compact roll of sterile cotton gauze that has excellent fluid absorption.	
8		Sterile non-woven gauze	7.5cm × 7.5cm-2pcs	10	рс	Appropriate for wound dressing, cleaning, prepping, or packing.	And a set of the set o
9		Wooden tongue depressor	16cm	5	ea	Used to depress the tongue during examination; during emergencies, insert it into the molar area, gently rotate, and separate the upper and lower teeth.	
10		Mouth opener	/	2	ea	Used to slowly rotate in order to open the oral cavity during emergencies.	



NO	MODULE	NAME	SPECFICATIONS	QTY	PIECE	FUNCTION	PICTURE
11		Tongue forceps	16cm	1	bu	Used during emergencies to place the tongue forceps at the base of the patient's tongue, gently clamp the tongue, and gently pull it out to ensure airway patency.	8
12		Guedel airway	3#	3	рс	Maintain or open the airway by stopping the tongue from covering the epiglottis.(Use an oropharyngeal airway only if the patient is unconscious or minimally responsive).	
13		Nasopharyngeal airway (including lubricant)	7#	3	рс	Assist with patient oxygenation and ventilation.	
14		Simple respirator	L	1	set	Used to provide lung ventilation during emergency artificial respiration for patients experiencing sudden respiratory distress or respiratory failure.	
15		Disposable laryngeal mask (including lubricant)	4#	1	рс	Used to maintain an open airway, prevent vomit from entering the respiratory tract, and facilitate surgery.	



NO	MODULE	NAME	SPECFICATIONS	QTY	PIECE	FUNCTION	PICTURE
16		Thoracentesis needle	14G	6	рс	Puncturing the pleural cavity with a puncture needle to drain accumulated air or fluid from within.	
17		Chest seal tape	15cm × 15cm	2	рс	designed to treat open chesttrauma which may lead to traumatictension pneumothorax under mostcircumstances.	
18		Bandage triangular	96cm × 96cm × 136cm	4	ea	One of the most versatile first aid dressings that come with two pins.	Bandage Pringular Transverter Transverter C
19		Pelvic Sling	87cm × 20cm	1	рс	In the case of a patient with pelvic fracture resulting in pelvic instability, a pelvic sling is used as a temporary stabilizing device to prevent further serious injury to the pelvis or severe bleeding.	
20		Eye shield	Hard eye patch + sterile eye sticker	1	ea	Used for non-chronic wound care, providing a microenvironment to facilitate wound healing.	



NO	MODULE	NAME	SPECFICATIONS	QTY	PIECE	FUNCTION	PICTURE
21		Elastic bandage	7.5cm × 450cm	6	rl	Delivers consistent pressure, making it ideal for situations that require higher compression.	
22		Self adhesive elastic bandage	10cm × 450cm	2	rl	Used for wound dressings or providing compression to limbs, serving the purpose of dressing and fixation.	
23		Medical external fixation splint	11cm × 46cm	2	рс	Can be used easily to splint any part of the body. Just mold it to shape, wrap or strap it into place.	
24		Medical Tape	2.5cm × 500cm	2	rl	Single use tape for securing first aid dressings, bandages and splints.	
25		Intravenous infusion set	three-port valve 1pc	1	set	Using the principles of atmospheric pressure and hydrostatic pressure to deliver large volumes of sterile liquids, electrolytes, and medications into the body via veins.	



NO	MODULE	NAME	SPECFICATIONS	QTY	PIECE	FUNCTION	PICTURE
26		Intravenous infusion set	Single intravenous indwelling needle 18G 1pc			Using the principles of atmospheric pressure and hydrostatic pressure to deliver large volumes of sterile liquids, electrolytes, and medications into the body via veins.	
27		Intravenous infusion set	single-use syringe 10ml 1pc			Using the principles of atmospheric pressure and hydrostatic pressure to deliver large volumes of sterile liquids, electrolytes, and medications into the body via veins.	
28		Intravenous infusion set	sterile non-woven gauze 7.5cm × 7.5cm,2pcs;1pc			Using the principles of atmospheric pressure and hydrostatic pressure to deliver large volumes of sterile liquids, electrolytes, and medications into the body via veins.	Constant of the second
29		Intravenous infusion set	Single-use flow setting fine infusion set 1pc			Using the principles of atmospheric pressure and hydrostatic pressure to deliver large volumes of sterile liquids, electrolytes, and medications into the body via veins.	
30		Intravenous infusion set	TPE tourniquet 2.5cm × 45cm 1pc			Using the principles of atmospheric pressure and hydrostatic pressure to deliver large volumes of sterile liquids, electrolytes, and medications into the body via veins.	0



NO	MODULE	NAME	SPECFICATIONS	QTY	PIECE	FUNCTION	PICTURE
31		Intravenous infusion set	Lien needle stick 1pc			Using the principles of atmospheric pressure and hydrostatic pressure to deliver large volumes of sterile liquids, electrolytes, and medications into the body via veins.	
32		Intravenous infusion set	Povidone-lodine Rrep Pad 6cm × 3cm, 2pcs			Using the principles of atmospheric pressure and hydrostatic pressure to deliver large volumes of sterile liquids, electrolytes, and medications into the body via veins.	Antingen Povidone-todino Prop Pad CE
33		Disposable syringes	2ml	8	рс	This product is used for clinical suction of liquids or for immediate injection after liquid infusion.	ľ
34			5ml	8	рс	This product is used for clinical suction of liquids or for immediate injection after liquid infusion.	
35			10ml	5	рс	This product is used for clinical suction of liquids or for immediate injection after liquid infusion.	and the second sec



NO	MODULE	NAME	SPECFICATIONS	QTY	PIECE	FUNCTION	PICTURE
36		Povidone-Iodine Rrep Pad	6cm × 3cm	10	рс	The wipes are alcohol free and sting free, which is suitable for cleaning up the wounds before applying dressings or other treatments.	Antenets Poylidone-Jodine Prep Pad CE 15 Transmission and Frequencies
37		Gauze bandage scissors	15cm	1	bu	Bandage scissor makes bandage removal very easy, smooth, and quick.	
38		Biohazard Bag	61cm × 81cm	2	bu	Used for handling biohazardous waste.	
39		Flashlight	134mm x 15mm x 15mm	1	ea	Perfect for inspection work and other emergency uses.	
40		Emergency blanket	210cm × 160cm	2	рс	It is widely used in emergencies and disasters to reduce heatloss.	A construction of the second s



NO	MODULE	NAME	SPECFICATIONS	QTY	PIECE	FUNCTION	PICTURE
41		Disposable single use sterilized latex examine glovel	Μ	5	ра	Reduce infection risk when dressing wounds or dealing withany body fluid or waste.	
42		Grinding wheel	Φ20mm	2	рс	Used for opening medical ampoules.	$\overline{\mathbf{\cdot}}$
43		Injury classification card	4 sheets/pack	3	рс	Rapid injury identification	1 1 winds were 1 winds 1 win
44		Injury ticket	50 pcs	1	рс	Recording casualty case documents.	INJURY TICKET
45		Marking pen	Double-ended purple	2	рс	Versatile dual tip design is excellent for broad writing on large spaces and ultra-fine precision marking on small surfaces.	



NO	MODULE	NAME	SPECFICATIONS	QTY	PIECE	FUNCTION	PICTURE
46		Medicine box	7.5cm × 6.3cm × 2.5cm	1	ea	Used for storing pills or tablets in a small container.	
47		Soft stretcher	180cm × 70cm	1	ea	Suitable for emergency rescue and short-distance transport.	
48		Bag	330x180x510mm	1	ea		

CONTACT

For inquiries regarding our product offerings, pricing structures, or contractual terms, we encourage you to reach out to our team.

European Office:

E-mail: contact@resqsupply.com

Phone: +45 911 19 911



Jaegersborg Allé 93 2820 Gentofte Denmark

CVR/VAT: 44646021

